


## Total Oxylipins



### Deliverables

This pre-validated panel includes more than 126 oxylipins and offers **fully quantitative results** (expressed in pmol). In addition, the amount of each sample used is provided (weight or volume). This analysis package allows to measure total amount of oxylipins including both esterified and free form.



### Method

Oxylipins are analysed by **LC/ESI-MS/MS** with an Agilent 1290 HPLC system combined with an Agilent 6495 Triplequad mass spectrometer (Agilent Technologies, Santa Clara, USA). Samples are extracted with organic solvents and alkaline hydrolysis is performed. For lipid quantification, deuterated internal standards are added prior to sample extraction. Crude extracts are further purified by solid phase extraction (SPE). Clean extracts are loaded on the LC system for LC/ESI-MS/MS analysis. Mass spectrometry analysis is performed with multiple reaction monitoring (MRM) in negative ion mode with at least two mass transitions for each compound. Data analysis is done using Mass Hunter software (Agilent Technologies).



### Sample submission guidelines

SAMPLE TYPE	SAMPLE AMOUNT	SHIPMENT CONDITIONS
<b>Whole blood, plasma, serum</b>	500 µl	<i>Frozen, on dry ice</i>
<b>Tissue</b>	40 mg	<i>Frozen pieces of tissues, on dry ice</i>
<b>Feces</b>	20 mg	<i>Frozen, on dry ice</i>
<b>Cells</b>	3 × 10 <sup>6</sup> cells or 500 µg of total protein	<i>Frozen cell pellet, on dry ice</i>
<b>Cell supernatant, media, liquids, etc.</b>	0.5–20 ml	<i>Frozen, on dry ice</i>

**Minimal sample number: 10**

**Standard turnaround time: 6 weeks**

(for projects with up to 100 samples)

**Expedited delivery** is possible for an additional fee.

### Applications

- > Nutrition (omega-3/6)
- > Inflammation
- > CYP/LOX pathway
- > Lipid hydroperoxide generation & turnover
- > Lipid oxidation

## Covered analytes

ANALYTE	FATTY ACID	CATEGORY	SUBCLASS	LLOQ* [pg]
10,11-DiHDPA	DHA	Docosanoids		5.3
10,11-EDP	DHA	Docosanoids		2.7
10-HDHA	DHA	Docosanoids		2.9
11,12-DHET	ARA	Eicosanoids		2.1
11,12-DiHETE	EPA	Eicosanoids		2.8
11,12-EEQ	EPA	Eicosanoids		6.5
11,12-EET	ARA	Eicosanoids		2.5
11-HDHA	DHA	Docosanoids		3.5
11-HEPE	EPA	Eicosanoids		52.7
11-HETE	ARA	Eicosanoids		5.1
12,13-DiHOME	LA	Octadecanoids		1.2
12,13-EpOME	LA	Octadecanoids		5.1
12-HEPE	EPA	Eicosanoids		4.6
12-HETE	ARA	Eicosanoids		6.7
12-HETrE	DGLA	Eicosanoids		4.9
12-HHTrE	ARA	Eicosanoids		2.3
12-HpEPE	EPA	Eicosanoids		2.3
12-HpETE	ARA	Eicosanoids		2.9
12-oxo-ETE	ARA	Eicosanoids		46.4
12-oxo-LTB4	ARA	Eicosanoids	Leukotriene	12.1
13,14-DiHDPA	DHA	Docosanoids		4
13,14-EDP	DHA	Docosanoids		3.5
13-gamma-HOTrE	GLA	Octadecanoids		3.2
13-HDHA	DHA	Docosanoids		15.6
13-HODE	LA	Octadecanoids		3
13-HOTrE	ALA	Octadecanoids		3.2
13-HpODE	LA	Octadecanoids		5.3
13-HpOTrE	ALA	Octadecanoids		4.9
13-HpOTrE-gamma	GLA	Octadecanoids		24.3
13-oxo-ODE	LA	Octadecanoids		1.9
14,15-DHET	ARA	Eicosanoids		2.6
14,15-DiHETE	EPA	Eicosanoids		3.9
14,15-EEQ	EPA	Eicosanoids		10.3
14,15-EET	ARA	Eicosanoids		2.6
14,15-EpEDE	DGLA	Eicosanoids		3.3
14-HDHA	DHA	Docosanoids		2.8
15(R)epi-LXA4	ARA	Eicosanoids	Lipoxin	6.2
15-HEPE	EPA	Eicosanoids		2.9

ANALYTE	FATTY ACID	CATEGORY	SUBCLASS	LLOQ [pg]
15-HETE	ARA	Eicosanoids		5.1
15-HETrE	DGLA	Eicosanoids		5.1
15-HpEPE	EPA	Eicosanoids		4.6
15-HpETE	ARA	Eicosanoids		4.6
15-oxo-ETE	ARA	Eicosanoids		6.9
16,17-DiHDPA	DHA	Docosanoids		4.7
16,17-EDP	DHA	Docosanoids		34.9
16-HDHA	DHA	Docosanoids		2.9
16-HETE	ARA	Eicosanoids		2.8
17(R)RvD1	DHA	Docosanoids	Resolvin	6.3
17(R,S)-RvD4	DHA	Docosanoids	Resolvin	5.9
17,18-DiHETE	EPA	Eicosanoids		6
17,18-EEQ	EPA	Eicosanoids		7.8
17-HDHA	DHA	Docosanoids		9.4
17-HETE	ARA	Eicosanoids		2.2
17-HpDHA	DHA	Docosanoids		3.3
18-COOH-dinor-LTB4	ARA	Eicosanoids	Leukotriene	8.9
18-HEPE	EPA	Eicosanoids		5.2
18-HETE	ARA	Eicosanoids		2.8
19,20-DiHDPA	DHA	Docosanoids		6.3
19,20-EDP	DHA	Docosanoids		5.8
19-HEPE	EPA	Eicosanoids		9
19-HETE	ARA	Eicosanoids		7.3
20-COOH-ARA	ARA	Eicosanoids		26.2
20-COOH-LTB4	ARA	Eicosanoids	Leukotriene	8.9
20-HDHA	DHA	Docosanoids		4.4
20-HEPE	EPA	Eicosanoids		31.4
20-HETE	ARA	Eicosanoids		25.4
22-HDHA	DHA	Docosanoids		8.7
4-HDHA	DHA	Docosanoids		3.3
4-oxo-DHA	DHA	Docosanoids		32.3
5(S),15(S)-DiHETE	ARA	Eicosanoids		12.8
5(S),12(S)-DiHETE	EPA	Eicosanoids		4
5,6-DHET	ARA	Eicosanoids		2
5,6-DiHETE	EPA	Eicosanoids		3.2
5,6-EEQ	EPA	Eicosanoids		250
5,6-EET	ARA	Eicosanoids		19.6
5-HEPE	EPA	Eicosanoids		4.1
5-HETE	ARA	Eicosanoids		4.1
5-HETrE	DGLA	Eicosanoids		2.1

ANALYTE	FATTY ACID	CATEGORY	SUBCLASS	LLOQ [pg]
<b>5-HpEPE</b>	EPA	<i>Eicosanoids</i>		5
<b>5-HpETE</b>	ARA	<i>Eicosanoids</i>		5
<b>5-oxo-ETE</b>	ARA	<i>Eicosanoids</i>		9.5
<b>6(S)-LXA4</b>	ARA	<i>Eicosanoids</i>	<i>Lipoxin</i>	6.3
<b>6-trans-epi-LTB4</b>	ARA	<i>Eicosanoids</i>	<i>Leukotriene</i>	4.2
<b>6-trans-LTB4</b>	ARA	<i>Eicosanoids</i>	<i>Leukotriene</i>	4.2
<b>7,8-DiHDPA</b>	DHA	<i>Docosanoids</i>		2.8
<b>7,8-EDP</b>	DHA	<i>Docosanoids</i>		3.9
<b>7-epi-Mar1</b>	DHA	<i>Docosanoids</i>	<i>Maresin</i>	5.8
<b>7-HDHA</b>	DHA	<i>Docosanoids</i>		5.5
<b>8(S),15(S)-DiHETE</b>	ARA	<i>Eicosanoids</i>		17.7
<b>8,9-DHET</b>	ARA	<i>Eicosanoids</i>		5.2
<b>8,9-DiHETE</b>	EPA	<i>Eicosanoids</i>		3.6
<b>8,9-EEQ</b>	EPA	<i>Eicosanoids</i>		10.4
<b>8,9-EET</b>	ARA	<i>Eicosanoids</i>		57.2
<b>8-HDHA</b>	DHA	<i>Docosanoids</i>		10.7
<b>8-HEPE</b>	EPA	<i>Eicosanoids</i>		4.3
<b>8-HETE</b>	ARA	<i>Eicosanoids</i>		49.2
<b>8-HETrE</b>	DGLA	<i>Eicosanoids</i>		9.2
<b>9,10-DiHOME</b>	LA	<i>Octadecanoids</i>		1.1
<b>9,10-EpOME</b>	LA	<i>Octadecanoids</i>		5.4
<b>9-HEPE</b>	EPA	<i>Eicosanoids</i>		9.1
<b>9-HETE</b>	ARA	<i>Eicosanoids</i>		46.9
<b>9-HODE</b>	LA	<i>Octadecanoids</i>		5.9
<b>9-HOTrE</b>	ALA	<i>Octadecanoids</i>		8.6
<b>9-HpODE</b>	LA	<i>Octadecanoids</i>		6.4
<b>9-HpOTrE</b>	ALA	<i>Octadecanoids</i>		6.4
<b>9-oxo-ODE</b>	LA	<i>Octadecanoids</i>		11.1
<b>9-oxo-OTrE</b>	ALA	<i>Octadecanoids</i>		5.1
<b>LTB3</b>	ETA	<i>Eicosanoids</i>	<i>Leukotriene</i>	7.1
<b>LTB4</b>	ARA	<i>Eicosanoids</i>	<i>Leukotriene</i>	2.8
<b>LTB5</b>	EPA	<i>Eicosanoids</i>	<i>Leukotriene</i>	3
<b>LXA4</b>	ARA	<i>Eicosanoids</i>	<i>Lipoxin</i>	5.9
<b>LXA5</b>	EPA	<i>Eicosanoids</i>	<i>Lipoxin</i>	5.2
<b>LXB4</b>	ARA	<i>Eicosanoids</i>	<i>Lipoxin</i>	6.3
<b>Mar 1</b>	DHA	<i>Docosanoids</i>	<i>Maresin</i>	5.8
<b>Mar 2</b>	DHA	<i>Docosanoids</i>	<i>Maresin</i>	5.8
<b>NPD1</b>	DHA	<i>Docosanoids</i>		3.6
<b>NPDx</b>	DHA	<i>Docosanoids</i>		3.6

ANALYTE	FATTY ACID	CATEGORY	SUBCLASS	LLOQ [pg]
<b>RvD1</b>	DHA	<i>Docosanoids</i>	<i>Resolvin</i>	6.5
<b>RvD2</b>	DHA	<i>Docosanoids</i>	<i>Resolvin</i>	5.6
<b>RvD3</b>	DHA	<i>Docosanoids</i>	<i>Resolvin</i>	5.3
<b>RvD4</b>	DHA	<i>Docosanoids</i>	<i>Resolvin</i>	4
<b>RvD5</b>	DHA	<i>Docosanoids</i>	<i>Resolvin</i>	3.9
<b>RvE1</b>	EPA	<i>Eicosanoids</i>	<i>Resolvin</i>	3.7
<b>RvE2</b>	EPA	<i>Eicosanoids</i>	<i>Resolvin</i>	2
<b>RvE4</b>	EPA	<i>Eicosanoids</i>	<i>Resolvin</i>	2
<b>tetranor-12-HETE</b>	ARA	<i>Eicosanoids</i>		2.2

\* lower limit of quantification

**For enquiries, please contact:**

Lipotype GmbH | Tatzberg 47 | 01307 Dresden | Germany  
 Lipotype Inc. | 1 Broadway, 14th floor | Cambridge, MA 02142 | United States  
 sales@lipotype.com



[www.lipotype.com](http://www.lipotype.com)

